

9F, No. 108, Jhouzih St.,Taipei, Taiwan Tel: + 886-2-8751-1888 Fax: + 886-2-6602-1218 E-mail: sales@abnova.com

Datasheet

COVID-19 Viral Antigen Rapid Test (CE IVD)

Catalog Number: DC0302

Regulatory Status: For research use only (RUO)

Product Description: COVID-19 Viral Antigen Rapid Test provides a fast, convenient, and qualitative testing of SARS nucleocapsid antigens from nasal swabs in individuals who are suspected of COVID-19 (SARS-CoV-2) within first 12 days of symptom onset.

Applications: LFIA

(See our web site product page for detailed applications information)

Protocols: See our web site at http://www.abnova.com/support/protocols.asp or product page for detailed protocols

Supplied Product: 25 Test Kit:

COVID-19 Antigen Testing Strips (25)

- Colloidal gold-labeled COVID-19 nucleocapsid murine monoclonal antibody
- COVID-19 nucleocapsid human monoclonal antibody
- Goat anti-murine polyclonal antibody

Extraction Buffers (25)

- Extraction reagent

Droppers (25)

- 100 uL Dropper for drawing sample mixture

Positive Control Swab (1)

- Non-infectious, recombinant viral nucleocapsid protein with less than 0.1% sodium azide

Negative Control Swab (1)

- Buffer with less than 0.1% sodium azide

Assay Documentation (1)

- Instructions for Use

Storage Instruction: Store product at 15°C-30°C (59°F-86°F) and avoid direct exposure to sunlight. Do not freeze or store the product outside the temperature range described above. Do not open until ready to use. Do not use it after the expiration date.